BODYJET

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Heart rate and oxygen saturation are measured once at the start of the program.

For accurate results, wear the sensor on your finger and remain still for at least one minute.

If the initial measurement is outside the normal range (Heart Rate: 60~100 bpm, Oxygen Saturation:

95~100%), a second measurement will be taken at the end of the program

to check whether your body functions have returned to the normal range after care.

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Cure - A care program designed to help improve chronic conditions such as colds, high blood pressure, and diabetes.

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Healing - Relieves mental fatigue such as stress and mild depression, while easing minor muscle aches.

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Sport - Promotes muscle development and removes lactic acid for effective post-workout recovery.

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Detox - Enhances lymphatic circulation and targets visceral fat to support full-body detox and weight loss.

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Silver - For seniors and patients: improves scalp circulation, supports mental well-being, and strengthens weakened muscles.

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Kids - For children aged 6 and above: enhances concentration and immunity while reducing electromagnetic stress.

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Free - A customizable mode allowing flexible, user-defined care settings.

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
Our company has preset default values for all programs.

These settings are based on clinical data collected from approximately 20,000 customers using the previous model over the past six years.

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Connect the steam hose to the main unit and secure it with the butterfly clamp.

※ Ensure the hose does not sag downward to prevent unnecessary malfunctions.

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Connect the steam generator's power cable and LED cable.

※ Check the number of pins on the ports and connect the 4-pin and 2-pin ports to the correct positions.

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Connect the drain nozzle to the external water container using the drain hose.

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Place the BodyJet cushion inside the main unit.

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Fill the steam generator with 2.5 liters of purified water and close the lid.

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Cover the steam outlet with a towel, select your desired program, and preheat for the set time.

Preheating time by program: Cure 25min / Healing 15min / Sport 20min / Detox 15min / Silver 25min / Kids 25min / Free 15-20min

* These times are standard; you may start once steam begins to flow during preheating.

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Enter the main unit after preheating.

※ Hold the handrails carefully to prevent injury.

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
Adjust the seat height according to your height.

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Wrap a towel around your neck to prepare for the session.

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
After the session, fully open the door, wipe off sweat while seated, and exit slowly.

※ Hold the handrails carefully to prevent injury.

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Height-Adjusting Footrest: Helps shorter adults or children use the unit more securely.

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Backrest Cover: For users with a slim body type, using the backrest allows for a more comfortable posture during the session.

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Press the [Settings] button at the bottom left to access the settings screen.

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
Enter the [PROGRAM SETTING] menu and select the program you want to adjust.

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TIME – Set the operation time of the steam generator. *Up to 45 minutes.

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LEVEL – Adjust the steam temperature within levels 8–16.

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RECLINE TIME – Set the time when the reclining function starts during the program.

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UPRIGHT TIME – Set the time when the reclining chair returns to its original position during the program.

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OK (Save) – Save your changes. ※Unsaved changes will not be applied.

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Level 12: approx. 32°C-37°C

Level 10: approx. 30°C-31°C

Level 8: approx. 28°C-29°C (resting level, no steam output)

“The steam temperature may vary depending on your body temperature and external environmental conditions.”

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
*Steam Focus by Level

Levels 14-16: Focuses on upper body

Levels 11-12: Focuses on abdomen and genital area

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

Program Setting Examples

1. If you set **TIME** to 30 minutes and **LEVEL** to 16, the device will operate for a total of 30 minutes and heat up to Level 16 (approximately 44°C-50°C).
2. If you then set it to 20 minutes at Level 14, this configuration will take effect **after the first 10 minutes** of the initial setting (30 minutes · Level 16).
In other words, the device will preheat at the highest temperature for the first 10 minutes, and when the customer enters at the 20-minute point, it will start at **Level 14 (approximately 38°C-43°C)**.
3. If you set it to 10 minutes at Level 12, the device will remain at Level 14 until there are 10 minutes remaining, and then switch to **Level 12 (approximately 32°C-37°C) at the 10-minute mark**.
4. If you then set it to 5 minutes at Level 8, the device will stay at Level 12 until there are 5 minutes remaining, and then lower the temperature to **Level 8 (approximately 28°C-29°C, resting level)**.
5. Finally, if you set it to **Level 16 for the last 2 minutes**, the temperature will rise again to the highest level just before the program ends.

Reclining Function Settings

6. The reclining start time must be set based on the moment the customer enters.
For example, if the customer enters when 20 minutes remain, you should not set the start time at “20 minutes.” Instead, adjust it to the desired time after the customer has entered.
*Currently, the default setting is 5 minutes after entry (i.e., at the 15-minute point).
7. The reclining return time is set based on the remaining time. For example, if the return time is set to 5 minutes, the chair will automatically return to its original position when 5 minutes remain.

If both the reclining start time and return time are set to “0,” the device will not move automatically.

In this case, you can use the manual reclining button during the session to adjust the backrest at any desired timing.

Free Mode Guide

- By entering the Free Mode program settings in the same way, users can customize and configure their own care programs.
- All programs other than Free Mode are pre-programmed in the device with optimal settings based on six years of clinical data (with approximately 20,000 participants).

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
Remove the inner pot from the steamer, pour out the water, and wipe the bottom of the pot with a dry towel.

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A/S Inquiry and Customer Service Center

(562) 529-5343

<http://www.umineusa.com>